



**Specifications:**

Dimensions (L x W x H):	22 x 7 x 14in. (500 x 190 x 360mm)
Weight:	44lbs.(20 kg)
Power requirements:	AC,90-240V, 50/60 Hz, 150 Watts
Laser model:	Ultraviolet laser
Laser power:	1 Watt
Working temperature:	59°F ~ 86°F (15°C ~ 30°C)
Work noise:	≤65dB
Connectivity:	Type.B
Cassette printing speed:	< 2.5 seconds per cassette (on-demand printing)
High resolution:	500 DPI
Number of hoppers :	1 hopper for cassettes with or without lids
Cassette capacity per hopper:	80 cassettes per hopper
Total cassette capacity:	80 cassettes
Output tray capacity:	8 cassettes

Guangdong Jinqun Medical tissue cassette laser printer is a must-have for pathology and lab environments. Precision construction forges jam-free printing while its low work noise ≤65dB design offers ribbon or cartridge-free printing, optimal for cost-efficient workflows.

Easily integrated into laboratory and hospital information systems, this laser cassette printer is perfectly compatible with cassettes with lid or without lid and Sakura's parofoam cassettes frames. Its structure possesses 1 cassette magazines holding up to 80 cassettes total. Choose from many printing options, including 1D as well as 2D barcoding, alphanumeric, graphic, and logo printing. On-demand mode printing, its high printing speed completes a cassette every 2.5 seconds.

**Features:**

- Compact size for saving footprint.
- Simple and powerful integrated into LIS system.
- Fast printing speed up to 2.5 seconds.
- Crisp and permanent printout resistant to processed solvents.
- Collector is flexible to face toward left or right.

Model No.	Description
EMR810	Laser Cassette Printer

